


Azithromycin

Experimental

Column:	SVEA C18 Opal, 5 µm 110 Å 4.6 x 250 mm
Instrument:	HPLC
Mobile phase:	Mobil A: 1.8g/L anhydrous disodium hydrogen phosphate adjusted till pH=8.9 by dilute phosphoric acid or dilute sodium hydroxide. Mobil B: methanol: acetonitrile (250:750 V/V). Gradient condition see following table.
Flow rate:	1.0 ml/min
Injection volume:	50µl
Sample concentration:	8 mg/ml
Column temperature:	60°C
Detector:	UV 210 nm
Sample:	Azithromycin including impurities



Peak order from left to right : M, Q, R, F, J, I, S, H, unknown, A, unknown, azithromycin, B

Gradient condition:

Time (mins)	Mobil phase A	Mobil phase B
0-25	50-45	50-55
25-30	45-40	55-60
30-80	40-25	60-75
80-81	25-50	75-50
81-93	50	50

Method: EU Pharmacopoeia 9.0, Azithromycin monograph. Data kindly provided by Yunbo Technologies.

